medi

M.4[®]s

knee brace

INSTRUCTIONS FOR USE / INSTRUCCIONES PARA USO

M.4[®]s

Dear Patient,

Your doctor has prescribed the functional knee brace M.4^s for you. Please read through these instructions carefully to help ensure that your brace is fitted correctly.

Indications

- Ruptures of cruciate ligaments with accompanying injuries.
- Isolated ruptures of the anterior cruciate ligament.
- Ruptures and/or instability of the collateral ligaments.
- After reconstruction of cruciate ligaments.
- For conservative treatment of knee joint instability.

Side-Effects / Contraindications

None are currently known.

Setting the flexion / extension limits

Extension limit wedge, silver:

0°, 10°, 20°, 30°, 45°

Flexion limit wedge, black:

0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 120°

Can be immobilized at: 0°, 10°, 20°, 30°, 45°

For immobilizing the orthosis, first insert the flexion wedge and then the corresponding extension wedge. The medial (inside) and the lateral (outside) wedge settings must always be the same.

1. 10° extension and 90° flexion wedges are pre-installed in your knee brace.
2. To change the wedges, simply remove the joint cover.
3. Replacement extension wedges (silver) can be found in the wedge set marked 'R-E' and 'L-E'. Slide the new wedge into the joint from the front, and secure it by tightening the screws. Insert the wedge so that the letter shown on the joint plate points upwards (R or L). Now insert the screw and tighten.
4. To replace the black flexion wedge, loosen the screw and pull the wedge out from the back of the orthosis. Select the desired wedge from the set marked 'R-F' and 'L-F'. Insert the wedge so that the letter shown on the joint plate points upwards (R or L). Now

insert the screw and tighten.

5. Bend and flex the brace 'joint' a couple of times and check that the wedges are positioned correctly. Finally, secure the joint cover.

Important notes

Changes to the extension and flexion limits should only be made with instructions from your doctor or orthotist. In order to prevent hyperextension, always insert extension wedges or 0° wedges. The wedges define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert the larger wedge to limit the joint movement to the desired excursion. In addition, make sure that the orthosis joints are correctly positioned: the middle of the joints should be level with the upper border of the patella. The joint must lie behind the midline of the leg.

Shaping the orthosis

When shaping the orthosis, please ensure that the joints stay parallel. This guarantees that the orthosis will function well and run easily and that there is no excessive wear.

Fastening the belt eyes

The belt eyes have been adequately fastened and are, to some extent, deliberately mobile. Please do not tighten the screws. If too much force is applied, the screws may break.

Fitting instructions

- Open all the straps. Sit on the edge of a chair, pull the brace up the affected leg and flex your knee to 45°. Position the joints so that the center of the hinges (extension wedges) are level with the top edge of the kneecap. Make sure that the brace is not twisted on your leg.
- First, fasten the lower front calf strap (strap #1), then the upper back thigh strap (#2). Keep your knee in a flexed position and push both joints towards the back of the leg.
- Now fasten the lower back thigh strap (#3). Make sure that the hinge joints are not positioned in front of the midline of the leg.
- Now fasten the upper front calf strap (#4). Then the upper back calf strap (#5), and finally the front thigh strap (#6).
- Bend your knee slightly a few times and refasten the straps again in the order described above to prevent the brace slipping.

What can I do if...

... the brace slips?

- Check that the brace is fitted properly. The center of the brace joints should be level with the upper border of the kneecap.
- Make sure that the back thigh strap (#5) is above the largest part of the calf muscle.
- Take the covers off the knee pads and wear the pads directly on your skin. Watch for skin reactions; skin irritation may occur in people with very sensitive skin.
- Do not wear the brace over clothes.

... the pads have worn down?

- Replacement pads can be obtained from mediUSA. Your orthotist will be pleased to help.

Care instructions

The inner pads and the condyle pads can be washed by hand with a mild soap/detergent. Leave to dry naturally. When used normally the brace is maintenance-free. However, after using the brace in water, the hinges can be lubricated using Teflon™ spray. Regularly check that the joint screws are secure, tighten if necessary.

Storage

Keep the brace in a cool, dry place and do not expose to direct sunlight.

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