medi

M2 Active OA Knee Brace

Instructions for use. Mode d'emploi. Instrucciones de uso.

medi. I feel better.

1

5

2

6

3

7

4

8

9

13

10

14

11

15

12

16

Brace

Intended purpose

The M2Active OA Knee Brace is a brace for load relief which works using a three-point principle. The product is exclusively to be used for the orthotic fitting of the knee and only on unbroken skin or if any wounds have been properly covered.

Indications

- Medial or lateral (unicondylar) osteoarthritis of the knee
- Following cartilage surgery and/or cartilage regeneration for unilateral post-operative load relief, where necessary
- For unicondylar stress reactions / stress fractures of the upper tibia

Contra-indications

Bilateral (bicondylar) osteoarthritis of the knee

Setting the restrictions for flexion and extension (can only be carried out by an orthopedic technician)

Extension restriction wedges, silver (stretch restriction): -5° , 0° , 10° , 20° , 30° , 45°

Flexion restriction wedges, black (bend restriction) 0° , 10° , 20° , 30° , 45° , 60° , 75° , 90° , 120°

Can be immobilised at: 0° , 10° , 20° , 30° , 45°

Adjusting the level of flexion and extension restriction

1. Your knee brace comes with 0° extension wedges pre-installed.
2. To change the wedges, remove the brace joint from the brace.
3. Loosen the screw on the inside section and remove the screw and wedge.

4. Place the appropriate wedge for the front and/or back in the side opening and screw the wedge in tight.
5. Move the brace joint a few times to ensure that it can move freely and the wedges are secure.
6. Mount the brace joint on the brace.

Only change the extension and flexion restrictions under the orders of the doctor treating you. To avoid hyperextension, extension wedges or 0° wedges must always be used. The wedges set the exact radius of movement for the brace. If soft tissue conditions pose issues, then it may be necessary to use a larger wedge, due to soft tissue compression, in order to limit movement of the knee to the desired radius.

Putting on the brace

Before use:

- Ensure that both load relief mechanisms are in the starting position (close to the knee joint) and that the switch is set to the white dot. (Fig. 1)
- Open all quick-release clasps on the front of the brace, by turning them slightly. (Fig. 2)
- Open the velcro fastening on the brace. (Fig. 3)

Putting the brace on:

- Bend your leg at approx. 70° and pull the brace onto the injured leg. (Fig. 4).
- Pull the brace up your leg until the fabric comes to below the kneecap and ensure that the kneecap is positioned in the center of the opening. The center of the brace joint between the two lateral joint splints is now sitting around the same level as the upper edge of the kneecap.
- Please make sure that the brace has been placed on the leg straight and is not twisted. (Fig. 5)

- Now close both of the velcro fastenings on the upper part of the brace. (Fig. 6)
- Fasten the quick-release clasp for strap 1 first and then fasten the quick-release clasp for strap 2 (Fig. 7).
- Now fasten straps 3 and 4. (Fig. 8)
- Push and/or pull the load relief mechanism up or down, moving along the increasing number scale. (Attention: To ensure the load relief mechanism works properly, the lock needs to be released and then refastened after setting). An audible click is heard for every change in level (Fig. 9 and 10) to show that the brace has engaged. The greatest level of load relief is attained at level 9.

The orthopaedic technician has set the brace in advance and has established a benchmark value with you. Depending on your level of pain, you can adjust the mechanism to fit your needs.

- Take a few steps and ensure that the brace is sitting correctly on the leg. If necessary, correct the position of the brace again.

Taking the brace off:

- Deactivate the load relief mechanism by releasing the switch and pushing the mechanism towards the knee joint. (Figs. 11 and 12)
- Then set the switch back to the white dot.
- Undo all of the quick-release clasps. (Fig. 13)
- Open the velcro fastening on the brace and pull it off your leg, moving downwards. (Fig. 14)
- Close the velcro fastening and all quick-release clasps to ensure that no damage is done to the brace. (Fig. 15)

Note for technicians:

- Set the load relief system (strap system) so that the patient attains a good level of pain relief. To do this, with the brace in its starting position (slider is set to position 1), pull the straps tightly. Then set up the three-point system. Level 7/8 is considered the guideline starting level. The patient can increase or reduce the level of pressure. Excess strap can be trimmed as necessary once final adjustment has been made (Fig. 16)
- Allow the patient to wear the brace for a few minutes to ensure that the straps are providing an appropriate and comfortable level of tension.

Care instructions

Velcro fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Air dry only.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a cool, dry place.

Material composition

Polyamide, PU, polyester, lycra, aluminum

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

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**Important notes**

This medical device is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist

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