medi

M2 Active OA Light Knee Brace

Instructions for use. Mode d'emploi. Instrucciones de uso.

medi. I feel better.

1

5

9

2

6

10

3

7

4

8

M2 Active OA Light Knee Brace

Intended purpose

The M2 Active OA Light Knee Brace is a brace for load relief which works using a three-point principle. The product is exclusively to be used for the orthotic fitting of the knee and only on intact skin or any wounds have been properly covered.

Indications

Medial or lateral (unicondylar) osteoarthritis of the knee

Contra-indications

Bilateral (bicondylar) osteoarthritis of the knee

Setting the restrictions for flexion and extension (can only be carried out by an orthopaedic technician)

Extension restriction wedges, silver (stretch restriction): -5°, 0°, 10°, 20°, 30°, 45°

Flexion restriction wedges, black (bend restriction) 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 120°

Can be immobilised at: 0°, 10°, 20°, 30°, 45°

Adjusting the level of flexion and extension restriction

1. Your knee brace comes with 0° extension wedges pre-installed.
2. To change the wedges, remove the brace joint from the brace.
3. Loosen the screw on the inside section and remove the screw and wedge.
4. Place the appropriate wedge for the front and/or back in the side opening and screw the wedge in tight.

5. Move the brace joint a few times to ensure that it can move freely and the wedges are secure.
6. Mount the brace joint on the brace.

Only change the extension and flexion restrictions under the orders of the doctor treating you. To avoid hyperextension, extension wedges or 0° wedges must always be used. The wedges set the exact radius of movement for the brace. If soft tissue conditions pose issues, then it may be necessary to use a larger wedge, due to soft tissue compression, in order to limit movement of the knee to the desired radius.

Putting on the brace

Before each use:

- Open the straps on the front and the velcro fastenings on the brace. (Fig. 1)

Putting the brace on:

- Bend your leg at approx. 70° and pull the brace onto the injured leg. (Fig. 2)
- Pull the brace up your leg until the fabric comes to below the kneecap and ensure that the kneecap is positioned in the center of the opening. The center of the brace joint between the two lateral joint splints is now sitting around the same level as the upper edge of the kneecap. Please make sure that the brace has been placed on the leg straight and is not twisted. (Fig. 3)
- Now close the velcro fastenings on the upper part of the brace. (Fig. 4)
- Once again, ensure that the brace is sitting correctly.

To adjust unloading pressure:

- Pull the lower strap (1) through the slot on the joint and set it to the setting that was defined by the technician. Fasten the velcro. Then set the upper strap (2) and fasten it.

- The straps at the back have already been set by the technician. (Figs. 5 and 6)
- Take a few steps and ensure that the brace is sitting correctly on the leg. If necessary, correct the position of the brace again. (Fig. 7)

Taking the brace off:

- Open both of the straps and stick the end of the straps back onto their respective straps, to prevent them from snagging clothing. (Fig. 8)
- Open the velcro fastening on the brace and pull it off your leg, moving downwards. Close the velcro fastening to prevent damage from occurring. (Fig. 9)

Note for technicians:

- Set the load relief system (strap system) so that the patient attains a good level of pain relief. The base setting should be set so that the patient can still increase or decrease pressure. Excess strap material can be trimmed away as needed after initial fitting. (Fig. 10)
- Allow the patient to wear the brace for a few minutes to ensure that the straps are providing an appropriate and comfortable level of tension.

Care instructions

Velcro fasteners should be closed for washing. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Air dry only
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the product in a cool, dry place.

Material composition

Polyamide, poliurethane, polyester, lycra, aluminum

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

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**Important notes**

This medical device is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist

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