medi GmbH & Co. KG
Medicusstraße 1
D-95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
 83 Fennell Street
 North Parramatta NSW 2151
 Australia
 T +61-2 9890 8696
 F +61-2 9890 8439
 sales@mediaustralia.com.au
 www.mediaustralia.com.au

medi Austria GmbH
 Adamgasse 16/7
 6020 Innsbruck
 Austria
 T +43 512 57 95 15
 F +43 512 57 95 15 45
 vertrieb@medi-austria.at
 www.medi-austria.at

medi Belgium bvba
 Posthoornstraat 13/1
 3582 Koersel
 Belgium
 T: + 32-11 24 25 60
 F: +32-11 24 25 64
 info@medibelgium.be
 www.medibelgium.be

medi Brasil
 Rua Cristovan de Vita, n.º 260, Un. 30,
 Centro Logístico Raposo Tavares,
 Vila Camargo
 Vargem Grande Paulista - SP
 CEP 06730-000
 Brasil
 T: +55-11-3500 8005
 sac@medibrasil.com
 www.medibrasil.com

MAXIS a.s.,
 medi group company
 Slezská 2127/13
 120 00 Prague 2
 Czech Republic
 T: +420 571 633 510
 F: +420 571 616 271
 info@maxis-medica.com
 www.maxis-medica.com

medi Danmark ApS
 Vejlegardsvej 59
 2665 Vallensbaek Strand
 Denmark
 T +45-70 25 56 10
 F +45-70 25 56 20
 kundeservice@sw.dk
 www.meditandmark.dk

medi Bayreuth Espana SL
 C/Canigo 2-6 bajos
 Hospitalet de Llobregat
 08901 Barcelona
 Spain
 T +34-932 60 04 00
 F +34-932 60 23 14
 medi@mediespana.com
 www.mediespana.com

medi France
 Z.I. Charles de Gaulle
 25, rue Henri Farman
 93297 Tremblay en France Cedex
 France
 T +33-1 48 61 76 10
 F +33-1 49 63 33 05
 infos@medi-france.com
 www.medi-france.com

medi Hungary Kft.
 Bokor u. 21.
 1037 Budapest
 Hungary
 T +36 1 371-0090
 F +36 1 371-0091
 info@medi.hu
 www.medi.hu

medi Nederland BV
 Heusing 5
 4817 ZB Breda
 The Netherlands
 T +31-76 57 22 555
 F +31-76 57 22 565
 info@medi.nl
 www.medi.nl

medi Polska Sp. z o.o.
 Zygmunta Starego 26
 44-100 Gliwice
 Poland
 T +48-32 230 60 21
 F +48-32 202 87 56
 info@medi-polska.pl
 www.medi-polska.pl

medi Bayreuth Unipessoal, Lda
 Rua do Centro Cultural, no. 43
 1700-106 Lisbon
 Portugal
 T +351-21 843 71 60
 F +351-21 847 08 33
 medi.portugal@medibayreuth.pt
 www.medi.pt

medi RUS LLC
 Business Center NEO GEO
 Butlerova Street 17
 117342 Moscow
 Russia
 T +7-495 374 04 56
 F +7-495 374 04 56
 info@medirus.ru
 www.medirus.ru

medi Trading (Shanghai) Co. Ltd.
 Rm 1806, Building 1 Sandhill Plaza
 No. 2290 Zuchongzhi Road
 Pudong New District
 200131 Shanghai
 People's Republic of China
 T: +86-21 61761988
 F: +86-21 61769607
 infocn@medi.cn

Medi Turk Ortopedi Medikal Ic ve Dis Tic. AS
 Kültür mah. çaldıran sok.14/3
 06420 Çankaya/Ankara
 Turkey
 T +90 312 435 20 26
 F +90 312 434 22 82
 info@medi-turk.com
 www.medi-turk.com

medi UK Ltd.
 Plough Lane
 Hereford HR4 OEL
 Great Britain
 T +44-1432 37 35 00
 F +44-1432 37 35 10
 enquiries@mediuk.co.uk
 www.mediuk.co.uk

medi Ukraine LLC
 Tankova Str. 8, office 35
 Business-center „Flora Park“
 Kiev 04112
 Ukraine
 T +380 44 591 11 63
 F +380 44 392 73 73
 info@medi.ua
 http://medi.ua

medi USA L.P.
 6481 Franz Warner Parkway
 Whitsett, N.C. 27377-3000
 USA
 T +1-336 4 49 44 40
 F +1-888 5 70 45 54
 info@mediusa.com
 www.mediusa.com

medi PT control

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Important notes

This medical device is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist.

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medi PT control

Intended purpose

medi PT control is a patellar luxation brace. The product is exclusively to be used for the orthotic fitting of the knee and only on unbroken skin or if any wounds have been properly covered.

Indications

- Conservative or operative treatment of patella luxation in connection with a slight to medium knee instability
- Lateral misalignment in connection with a slight to medium knee instability
- Following lateral release operations
- Retropatellar cartilage surgery

Contraindications

None known to date.

Setting the flexion / extension limits Only to be carried out by an orthopedic technician.

Extension limit wedge, EXT (limits extension): 0°, 10°, 20°, 30°, 45°

Flexion limit wedge, FLEX (limits flexion): 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

Can be immobilised at:

0°, 10°, 20°, 30°, 45°

For immobilising the orthosis first insert the extension wedge and then the corresponding flexion wedge.

The medial (inside) and the lateral (outside) wedge settings must always be the same.

1. Your knee brace is provided with pre-installed 0° extension wedges.
2. In order to change the wedges, push the joint cover from the back toward the outside (see profile).
3. To change the wedges, remove the

tweezers from the joint cover and use them to simply pull the wedges toward the front (extension) or toward the back (flexion) out of the joint. Each wedge is labelled with „FLEX“ or „EXT“ and with the respective degree setting. Another distinguishing feature is the different colour of the wedges. Extension wedges are silver and flexion wedges are black. Insert the new wedge from the front (extension) or from the back (flexion) into the joint so that the „UP“ marking on the joint and on the wedge correspond. A clicking sound confirms that the wedge has interlocked.

4. Move the brace joint and test the position of the wedges. Now put back the tweezers into the holder and attach the joint cover again by hooking the two catches of the joint cover onto the front of the joint and then pushing the cover onto the joint. A clicking sound confirms that the joint cover is fastened securely.

Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist. In order to prevent hyperextension, always insert extension wedges or 0° wedges. The wedges define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert the larger wedge to limit the joint movement to the desired excursion. In addition, make sure that the orthosis hinges are correctly positioned: the middles of the hinges should be level with the upper border of the patella. The hinge must lie behind the midline of the leg.

Fitting instructions

- Open all the straps so that they will not interfere with you putting on the brace.

- Now put your leg into a 30-degree angle flexion position and pull the brace up over your leg (Image 1),
- Make sure that the kneecap is positioned centrally in the space for the patella. Position the brace so that the middles of the hinges are level with the upper border of the kneecap and just behind the lateral centre line of the leg (Image 2).
- Now close the straps in the following order:
 - Lower calf strap in the back (Image 3).
 - Lower thigh strap in the back (Image 4).
 - Patellar straps so that the kneecap is guided from the outside to the inside
 - please ensure that the patellar air cushion is positioned to the side of the kneecap and that it does not push onto the kneecap from the front (Image 5 + 6).
 - Upper thigh strap in the back (Image 7).
 - Thigh strap in the front (Image 8).
 - Calf strap in the front (Image 9).
- If possible, walk a few steps in order to check that the brace is in the correct position. Please adjust the position as required (Image 10).

Care instructions

Please close the Velcro fasteners before washing and remove the ortheses. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep in a dry place at room temperature. Protect from heat, direct sunlight and moisture.

Material composition

Avional, polyamide, elasthane, PU

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

wishes you a quick return to full fitness.

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