

medi Rx Knee Brace Under Sleeve

The makers of mediven medical and CEP Sports Compression products, with over 60 years of medical knitting expertise, offer a new knee brace undersleeve.

- Durable, high-quality medical grade, latex-free fabric for user comfort.
- Can be used with all knee braces.
- Silicone anti-migration top band.
- Anatomically designed for greater wearing comfort.
- Sizes extra small through double extra-large to accommodate the majority of application needs.



medi Rx Knee Brace Under Sleeve Sizing Chart

Size	Inches	Cm	Item Number
XS	12 1/4" - 14 1/2"	31 - 37 cm	GE09001
SM	14 1/2" - 17 3/4"	37 - 45 cm	GE09002
MD	17 3/4" - 20"	45 - 51 cm	GE09003
LG	20" - 22 1/2"	51 - 57 cm	GE09004
XL	22 1/2" - 25 1/2"	57 - 65 cm	GE09005
XXL	25 1/2" - 29"	65 - 74 cm	GE09006



*Measurement:
Circumference 6" (15cm)
above mid-patella.

Knee Brace Versions Available from medi USA



Single Upright OA in
Varus or Valgus versions



Functional Knee Brace in
Standard and Compact sizes



Double Upright OA available in
Varus and Valgus versions with
Standard and Compact sizes sizes