

protect.footsupports Junior



heat moldable

sports
everyday



The Junior insole has been developed specifically to provide maximum support and improved stability for feet from children through adults.

- Benefits**
- Deep heel cup for excellent support and control
 - Multi-layered core offers guidance and support of the foot
 - Slip-resistant bottom cover ensures excellent fit in the shoe
 - Heat moldable

- Indications**
- flat foot
 - instability

Size	item no
Euro 19/20, Kids 4-5	PIA1819
Euro 21/22, Kids 5.5-6.5	PIA1821
Euro 23/24, Kids 7-8	PIA1823
Euro 25/26, Kids 8.5-9.5	PIA1825
Euro 27/28, Kids 10-11	PIA1827
Euro 29/30, Kids 11.5-12.5	PIA1829
Euro 31/32, Kids 13-1	PIA1831
Euro 33/34, Youth 1.5-3	PIA1833
Euro 35/36, Youth 3.5-4.5 Women 4-5.5	PIA1835
Euro 37/38, Youth 5-6.5 Women 6-7.5	PIA1837
Euro 39/40, Youth 7-8 Women 8-9.5, Men 5-6.5	PIA1839
Euro 41/42, Youth 8.5-9.5 Women 10-11.5, Men 7-8.5	PIA1841
Euro 43/44, Youth 10-11 Women 12-13.5, Men 9-10.5	PIA1843
Euro 45/46, Youth 11.5-12.5 Men 11-12.5	PIA1845



Soccer shoes



Sneakers



Comfortable shoes



Gym shoes

medi USA
6481 Franz Warner Parkway
Whitsett, NC 27377
T 800-633-6334
F 888-570-4554
orders@mediusa.com
b2b.mediusa.com
mediusa.com

protect
footsupports



medi. I feel better.

protect.footsupports
Plantar Fasciitis Pro/ Plantar Fasciitis

sports everyday ●●●○○



Plantar Fasciitis Pro

Size	item no
Euro 35, US Men -, US Women 4-4.5	PIA0635
Euro 36, US Men -, US Women 5-5.5	PIA0636
Euro 37, US Men -, US Women 6-6.5	PIA0637
Euro 38, US Men -, US Women 7-7.5	PIA0638
Euro 39, US Men -, US Women 8-8.5	PIA0639
Euro 40, US Men 6-6.5, US Women 9-9.5	PIA0640
Euro 41, US Men 7-7.5, US Women 10-10.5	PIA0641
Euro 42, US Men 8-8.5, US Women 11-11.5	PIA0642
Euro 43, US Men 9-9.5, US Women 12-12.5	PIA0643
Euro 44, US Men 10-10.5, US Women -	PIA0644
Euro 45, US Men 11-11.5, US Women -	PIA0645
Euro 46, US Men 12-12.5, US Women -	PIA0646
Euro 47, US Men 13-13.5, US Women -	PIA0647
Euro 48, US Men 14-14.5, US Women -	PIA0648

The Plantar Fasciitis insole relieves pressure in the heel area and along the plantar fascia to ease foot pain.

Benefits

- Targeted soft padding of the heel area and along the plantar fascia
- Multi-layered core for support of the longitudinal arch
- Relief of painful areas
- Integrated metatarsal pad supports the transverse arch (Pro only)

Indications

- heel spur
- plantar fasciitis
- heel pain

Plantar Fasciitis

Size	item no
Euro 35/36, US Men -, US Women 4-5.5	PIA0535
Euro 37/38, US Men -, US Women 6-7.5	PIA0537
Euro 39/40, US Men -, US Women 8-9.5	PIA0539
Euro 41/42, US Men 7-8.5, US Women 10-11.5	PIA0541
Euro 43/44, US Men 9-10.5, US Women 12-12.5	PIA0543
Euro 45/46, US Men 11-12.5, US Women -	PIA0545
Euro 47/48, US Men 13-14.5, US Women -	PIA0547



protect.footsupports Active

heat moldable

sports everyday ●●●○○



The Active insole cushions feet during athletic activities with a targeted footbed with strategic support.

Benefits

- Optimal support for the transverse arch
- Quick and simple adaptation to virtually all sports shoes
- Multi-layered core for excellent support
- Heat moldable

Indications

- flat foot
- splay foot
- foot pain

Order information

Size	item no
Euro 35/36, US Men - US Women 4-5.5	PIA0435
Euro 37/38, US Men -, US Women 6-7.5	PIA0437
Euro 39/40, US Men -, US Women 8-9.5	PIA0439
Euro 41/42, US Men 7-8.5, US Women 10-11.5	PIA0441
Euro 43/44, US Men 9-10.5, US Women 12-12.5	PIA0443
Euro 45/46, US Men 11-12.5, US Women - 13-14	PIA0445
Euro 47/48, US Men 13-14.5, US Women - 14-15	PIA0447



protect.footsupports Comfort

heat moldable

sports everyday ●●○○○



The Comfort insole has been specially developed for feet that require extra cushioning such as diabetic or arthritic feet.

Benefits

- Pre-smoothed foundation for an individual fit
- Antibacterial top cover for a pleasant foot climate
- Multi-layered construction for a soft cushioned footbed
- Light core for support of the longitudinal arch
- Heat moldable

Indications

- diabetes
- arthritis
- sensitive skin

Order information

Size	item no
Euro 35/36, US Men -, US Women 4-5.5	PIA0135
Euro 37/38, US Men -, US Women 6-7.5	PIA0137
Euro 39/40, US Men -, US Women 8-9.5	PIA0139
Euro 41/42, US Men 7-8.5, US Women 10-11.5	PIA0141
Euro 43/44, US Men 9-10.5, US Women 12-12.5	PIA0143
Euro 45/46, US Men 11-12.5, US Women -	PIA0145
Euro 47/48, US Men 13-14.5, US Women -	PIA0147

